

A randomized crossover trial on the effect of compression stockings on nausea and vomiting in early pregnancy



BACKGROUND

50-80% of all women experience nausea and vomiting in early pregnancy (NVP). NVP has a detrimental impact on women's quality of life (QoL), especially their ability to work.

Despite the fact that a number of therapeutic interventions exist to alleviate symptoms – ranging from pharmaceutical products to natural therapies – women remain reluctant to take medications during pregnancy, and a safe, effective and acceptable treatment for NVP is still to be found.

STUDY DESIGN

The aim of this study was to evaluate the effect of medical compression stockings (MCS) on symptoms of nausea, vomiting and dizziness between 8 and 16 weeks of pregnancy, and their impact on QoL factors.

PARTICIPANTS


Pregnant women >18 years and with mild to moderate NVP were eligible to participate. A total of 74 women were enrolled in the study. 58 women completed the study and were included in the analysis set.



GROUP 1: «MCS first» = 2 weeks with // 2 weeks without

- Eligible: n=33
- Dropout: n=4 (3 miscarriages, 1 questionnaire incomplete)
- Final evaluation: n=29

PATIENT RANDOMIZATION




GROUP 2: «MCS second» = 2 weeks without // 2 weeks with


- Eligible: n=41
- Dropout: n=12 (3 miscarriages; 9 did not return for controls)
- Final evaluation: n=29



RESULTS: Endpoint = change in the Nausea and Vomiting in Pregnancy Quality of Life (NVPQOL) score




Compliance rate



Would recommend wearing MCS during pregnancy

With MCS



Significantly less nausea and vomiting



Fewer physical & physiological limitations; Less leg pain & dizziness



CONCLUSION

Wearing compression stockings in early pregnancy improves not only nausea and vomiting-associated symptoms, but also QoL factors such as dizziness, leg pain and physical/physiological limitations.

AUTHORS' SUGGESTION

Compression stockings could be a non-pharmaceutical instrument to alleviate nausea and vomiting symptoms during early pregnancy.

¹Reference: A randomized crossover trial on the effect of compression stockings in nausea and vomiting in early pregnancy. Mendoza E, Amsler F. Int J Womens Health. 2017;9: 89-99.