Compression Bulletin

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LITERATURVERZEICHNIS

- Mosti G, Iabichella ML and Partsch H. Compression therapy in mixed ulcers increases venous output and arterial perfusion. *J Vasc Surg* 2012; 55: 122–128.
- Dennis M, Cranswick G and Deary A. High-length versus below-knee stockings for deep venous thrombosis prophylaxis after stroke: a randomized trial. *Ann Intern Med* 2010; 153: 553–562.
- Reich-Schupke S, Feldhaus F, Altmeyer P, et al. Efficacy and comfort of medical compression stockings with low and moderate pressure six weeks after vein surgery. *Phlebologie* 2014; 29: 358–366.
- Harrison MB, Vandenkerkhof EG, Hopman WM, et al. The Canadian Bandaging Trial: evidence-informed leg ulcer care and the effectiveness of two compression technologies. *BMC Nurs* 2011; 10: 20.
- Mizuno J and In-Nami H. Allergic contact dermatitis to synthetic rubber, neoprene in compression stockings [In Japanese]. *Masui* 2011; 60: 104–106.
- Valesky EM, Kaufmann R and Meissner M. Contact allergy to compression stockings: is this possible? *Phlebologie* 2014; 43: 140–143.
- Lattimer CR, Azzam M, Kalodiki E, et al. Compression stockings significantly improve hemodynamic performance in post-thrombotic syndrome irrespective of class or length. *J Vasc Surg* 2013; 58: 158–165.
- Boris M, Weindorf S and Lasinski BB. The risk of genital edema after external pump compression for lower limb lymphedema. *Lymphology* 1998; 31: 15–20.
- Hansson C, Faergemann J and Swanbeck G. Fungal infections occurring under bandages in leg ulcer patients. *Acta Derm Venereol* 1987; 67: 341–345.
- Carlotta AB, Kaya E, Turgut H, et al. Letter to the editor. Folliculitis associated with intermittent pneumatic compression. *Yonsei Med J* 2014; 55: 545–546.
- O'Meara S, Al-Kurdi D, Ologun Y, et al. Antibiotics and antiseptics for venous leg ulcers. *Cochrane Database Syst Rev* 2014; CD003557.
- Franks PJ, Barker J, Collier M, et al. Management of patients with venous leg ulcers: challenges and current best practice. *J Wound Care* 2016; 25: S1–S67.
- O'Meara S, Cullum N, Nelson EA, et al. Compression for venous leg ulcers. *Cochrane Database Syst Rev* 2012; 11: CD000265.
- Elder DM and Greer KE. Venous disease: how to heal and prevent chronic leg ulcers. *Geriatrics* 1995; 50: 30–36.
- Dawn G, Loney M, Zamiri M, et al. Erosive pustular dermatosis of the leg associated with compression bandaging and fungal infection. *Br J Dermatol* 2003; 148: 489–492.
- Thompson CB, Wiemken TL and Brown TS. Effect of postoperative dressing on excisions performed on the leg: a comparison between zinc oxide compression dressings versus standard wound care. *Dermatol Surg* 2017; 43: 1379–1384.
- Jaccard Y, Singer E, Degischer S, et al. Effect of silver-threads-containing compression stockings on the cutaneous microcirculation: a double-blind, randomized cross-over study. *Clin Hemorheol Microcirc* 2007; 36: 65–73.
- Basford JR. The Law of Laplace and its relevance to contemporary medicine and rehabilitation. *Arch Phys Med Rehabil* 2002; 83: 1165–1170.
- Crawford F, Welch K, Andras A, et al. Ankle brachial index for the diagnosis of lower limb peripheral arterial disease. *Cochrane Database Syst Rev* 2016; 9: CD010680.
- Heath DI, Kent SJ, Johns DL, et al. Arterial thrombosis associated with graduated pressure antiembolic stockings. *Br Med J (Clin Res Ed)* 1987; 295: 580.
- Merrett ND and Hanel KC. Ischaemic complications of graduated compression stockings in the treatment of deep venous thrombosis. *Postgrad Med J* 1993; 69: 232–234.
- Creton D. La chirurgie des varices. Complications cutanées dues à la compression par doubles collants postopératoires [In French]. *Phlebologie* 1998; 51: 363–364.
- Ong JC, Chan FC and McCann J. Pressure ulcers of the popliteal fossae caused by thromboembolic deterrent stockings (TEDS). *Ir J Med Sci* 2011; 180: 601–602.
- Callam MJ, Ruckley CV, Dale JJ, et al. Hazards of compression treatment of the leg: an estimate from Scottish surgeons. *Br Med J (Clin Res Ed)* 1987; 295: 1382.
- Dennis M, Sandercock PA, Reid J, et al. Effectiveness of thigh-length graduated compression stockings to reduce the risk of deep vein thrombosis after stroke (CLOTS trial 1): a multicentre, randomised controlled trial. *Lancet* 2009; 373: 1958–1965.
- Chan CL, Meyer FJ, Hay RJ, et al. Toe ulceration associated with compression bandaging: observational study. *BMJ* 2001; 323: 1099.
- McIlhone S, Ukra H, Karim A, et al. Soft tissue injury to the sole of the foot secondary to a retained AV impulse foot pump. *Foot Ankle Surg* 2012; 18: 216–217.
- Anand A. Complications associated with intermittent pneumatic compression devices. *Anesthesiology* 2000; 93: 1556–1557.
- Parra RO, Farber R and Feigl A. Pressure necrosis from intermittent-pneumatic-compression stockings. *N Engl J Med* 1989; 321: 1615.
- Werbel GB and Shybut GT. Acute compartment syndrome caused by a malfunctioning pneumatic compression boot. A case report. *J Bone Joint Surg Am* 1986; 68: 1445–1446.
- Usmani N, Baxter KF and Sheehan-Dare R. Partially reversible common peroneal nerve palsy secondary to compression with four-layer bandaging in a chronic case of venous leg ulceration. *Br J Dermatol* 2004; 150: 1224–1225.
- Fukuda H. Bilateral peroneal nerve palsy caused by intermittent pneumatic compression. *Intern Med* 2006; 45: 93–94.
- McGrory BJ and Burke DW. Peroneal nerve palsy following intermittent sequential pneumatic compression. *Orthopedics* 2000; 23: 1103–1105.
- Hirate H, Sobue K, Tsuda T, et al. Peripheral nerve injury caused by misuse of elastic stockings. *Anaesth Intensive Care* 2007; 35: 306–307.
- Kim JH, Kim WI, Kim JY, et al. Peroneal nerve palsy after compression stockings application. *Saudi J Anaesth* 2016; 10: 462–464.
- Güzelkükük, Skempes D and Kummerde W. Common peroneal nerve palsy caused by compression stockings after surgery. *Am J Phys Med Rehabil* 2014; 93: 609–611.
- O'Brien C and Eltigani T. Common peroneal nerve palsy as a possible sequelae of poorly fitting below-knee thromboembolic deterrent stockings (TEDS). *Ann Plastic Surg* 2006; 57: 356–357.
- Blätler W and Partsch H. Leg compression and ambulation is better than bed rest for the treatment of acute deep venous thrombosis. *Int Angiol* 2003; 22: 393–400.

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39. Partsch H and Blättler W. Compression and walking versus bed rest in the treatment of proximal deep vein thrombosis with low molecular weight heparin. *J Vasc Surg* 2000; 32: 861–869.
40. Partsch H. The role of leg compression in the treatment of deep vein thrombosis. *Phlebology* 2014; 29: 66–70.
41. Ten Cate-Hoek AJ, Amin EE, Bouman AC, et al. Individualised versus standard duration of elastic compression therapy for prevention of post-thrombotic syndrome (IDEAL DVT): a multicentre, randomised, single-blind, allocation-concealed, non-inferiority trial. *Lancet Haematol* 2018; 5: e25–e33.
42. Partsch H. Therapy of deep vein thrombosis with low molecular weight heparin, leg compression and immediate ambulation. *VASA* 2001; 30: 195–204.
43. Scurr JH, Smith PD and Machin S. Deep vein thrombosis in airline passengers – the incidence of deep vein thrombosis and the efficacy of elastic compression stockings. *Cardiovasc Surg* 2001; 9: 159–161.
44. Jünger M, Partsch H, Kahle B, et al. Phlebologischer Kompressionsverband (PKV) – Leitlinie der Deutschen Gesellschaft für Phlebologie. *Phlebologie* 2009; 38: 168–171.
45. AWMF online. S3-Leitlinie Prophylaxe der venösen Thromboembolie (VTE), detail/II/003-001.html (2010, accessed July 2019).
46. AWMF online. S2k Leitlinie Diagnostik und Therapie der Lymphödem AWMF Reg. -Nr. 058-001, ien/058-001_S2k_Diagnostik_und_Therapie_der_Lymphoedeme_2017-05.pdf (2017, accessed July 2019).
47. Andriessen A, Apelqvist J, Mosti G, et al. Compression therapy for venous leg ulcers: risk factors for adverse events and complications, contraindications – a review of present guidelines. *J Eur Acad Dermatol Venereol* 2017; 31: 1562–1568.
48. Mostbeck A and Partsch H. Influence of dihydroergotamine and leg compression on the blood volume in different regions of the body [in German]. *Med Klin* 1978; 73: 801–806.
49. Lattimer CR, Kalodiki E, Azzam M, et al. Haemodynamic performance of low strength below knee graduated elastic compression stockings in health, venous disease, and lymphoedema. *Eur J Vasc Endovasc Surg* 2016; 52: 105–112.
50. Azevedo PS, Polegato BF, Minicucci MF, et al. Cardiac remodeling: concepts, clinical impact, pathophysiological mechanisms and pharmacologic treatment. *Arq Bras Cardiol* 2016; 106: 62–69.
51. Bain RJ, Tan LB, Murray RG, et al. Central haemodynamic changes during lower body positive pressure in patients with congestive cardiac failure. *Cardiovasc Res* 1989; 23: 833–837.
52. Dereppe H, Hoylaerts M, Renard M, et al. The effects of pressotherapy on cardiovascular hemodynamics [In French]. *Rev Med Brux* 1989; 10: 185–186.
53. Todd J, Austwick T, Berridge D, et al. B-type natriuretic peptide in lymphedema. *Lymphology* 2011; 44: 29–34.
54. The Criteria Committee of the New York Heart Association. Nomenclature and criteria for diagnosis of diseases of the heart and great vessels. Boston, MA: Little, Brown & Co, 1994, pp.253–256.
55. Galm O, Jansen-Genzel W, von Helden J, et al. Plasma human atrial natriuretic peptide under compression therapy in patients with chronic venous insufficiency with or without cardiac insufficiency. *VASA* 1996; 25: 48–53.
56. Wilputte F, Renard M and Venner JP. Hemodynamic response to multilayered bandages dressed on a lower limb of patients with heart failure. *Eur J Lym* 2005; XV: 1–4.
57. Leduc O, Cresset V, Leluc C, et al. Impact of manual lymphatic drainage on hemodynamic parameters in patients with heart failure and lower limb edema. *Lymphology* 2011; 44: 13–20.
58. Gorelik O, Almozino-Saraïan D, Litvinov V, et al. Seating-induced postural hypotension is common in older patients with decompensated heart failure and may be prevented by lower limb compression bandaging. *Gerontology* 2009; 55: 138–144.
59. Bousuges A, Ayme K, Chaumet G, et al. Observational study of potential risk factors of immersion pulmonary edema in healthy divers: exercise intensity is the main contributor. *Sports Med Open* 2017; 3: 35.
60. Carter HH, Spence AL, Ainslie PN, et al. Differential impact of water immersion on arterial blood flow and shear stress in the carotid and brachial arteries of humans. *Physiol Rep* 2017; 5: pii: 5/10/e13285.
61. Bowering CK. Use of layered compression bandages in diabetic patients. Experience in patients with lower leg ulceration, peripheral edema, and features of venous and arterial disease. *Adv Wound Care* 1998; 11: 129–135.
62. Robertson BF, Thomson CH and Siddiqui H. Side effects of compression stockings: a case report. *Br J Gen Pract* 2014; 64: 316–317.
63. Aschwanden M, Labs KH, Engel H, et al. Acute deep vein thrombosis: early mobilization does not increase the frequency of pulmonary embolism. *Thromb Haemost* 2001; 85: 42–46.
64. Schellong SM, Schwarz T, Kropp J, et al. Bed rest in deep vein thrombosis and the incidence of scintigraphic pulmonary embolism. *Thromb Haemost* 1999; 82: 127–129.
65. Partsch H and Mostbeck A. Early diagnosis of deep venous thrombosis of the lower leg [In German]. *Acta Med Austriaca* 1979; 4: 159–160.
66. Verlato F, Zucchetta P, Prandoni P, et al. An unexpectedly high rate of pulmonary embolism in patients with superficial thrombophlebitis of the thigh. *J Vasc Surg* 1999; 30: 1113–1115.
67. Boehler K, Kittler H, Stolkovich S, et al. Therapeutic effect of compression stockings versus no compression on isolated superficial vein thrombosis of the legs: a randomized clinical trial. *Eur J Vasc Endovasc Surg* 2014; 48: 465–471.
68. Abu-Own A, Shami SK, Chittenden SJ, et al. Microangiopathy of the skin and the effect of leg compression in patients with chronic venous insufficiency. *J Vasc Surg* 1994; 19: 1074–1083.
69. Dini V. Compression in vasculitis. *Veins and Lymphatics* 2016; 5: 5981.
70. Roumen-Klappé EM, den Heijer M, van Rossum J, et al. Multilayer compression bandaging in the acute phase of deep-vein thrombosis has no effect on the development of the post-thrombotic syndrome. *J Thromb Thrombolysis* 2009; 27: 400–405.
71. Amin EE, Joore MA, Ten Cate H, et al. Clinical and economic impact of compression in the acute phase of deep vein thrombosis. *J Thromb Haemost* 2018; 16: 1555–1563.
72. Wu SC, Crews RT, Najafi B, et al. Safety and efficacy of mild compression (18–25mm Hg) therapy in patients with diabetes and lower extremity edema. *J Diabetes Sci Technol* 2012; 6: 641–647.
73. Wu SC, Crews RT, Skratsky M, et al. Control of lower extremity edema in patients with diabetes: double blind randomized controlled trial assessing the efficacy of mild compression diabetic socks. *Diabetes Res Clin Pract* 2017; 127: 35–43.
74. Bochmann RP, Seibel W, Haase E, et al. External compression increases forearm perfusion. *J Appl Physiol* (1985) 2005; 99: 2337–2344.
75. Mayrovitz HN and Macdonald JM. Medical compression: effects on pulsatile leg blood flow. *Int Angiol* 2010; 29: 436–441.
76. Fromy B, Legrand MS, Abraham P, et al. Effects of positive pressure on both femoral venous and arterial blood velocities and the cutaneous microcirculation of the forefoot. *Cardiovasc Res* 1997; 36: 372–376.
77. Humphreys ML, Stewart AH, Gohel MS, et al. Management of mixed arterial and venous leg ulcers. *Br J Surg* 2007; 94: 1104–1107.
78. Alizadeh-Ghavidel A, Ramezannejad P, Mirmesdagh Y, et al. Prevention of edema after coronary artery bypass graft surgery by compression stockings. *Res Cardiovasc Med* 2014; 3: e17463.
79. Belczak CE, Tyska AL, Godoy JM, et al. Clinical complications of limb undergone harvesting of great saphenous vein for coronary artery bypass grafting using bridge technique. *Rev Bras Cir Cardiovasc* 2009; 24: 68–72.
80. Khoshgofar Z, Ayat Esfahani F, Marzban M, et al. Comparison of compression stocking with elastic bandage in reducing postoperative edema in coronary artery bypass graft patient. *J Vasc Nurs* 2009; 27: 103–106.
81. Maleti O. Compression after vein harvesting for coronary bypass. *Veins and Lymphatics* 2016; 5: 5989.
82. te Slaa A, Dolmans DE, Ho GH, et al. Prospective randomized controlled trial to analyze the effects of intermittent pneumatic compression on edema following autologous femoropopliteal bypass surgery. *World J Surg* 2011; 35: 446–454.
83. Olszewski WL. Episodic dermatolymphangioadenitis (DLA) in patients with lymphedema of the lower extremities before and after administration of benzathine penicillin: a preliminary study. *Lymphology* 1996; 29: 126–131.
84. Maccio A. Compression in dermatolymphangioadenitis. *Veins and Lymphatics* 2016; 5: 5982.
85. Stalbow J. Preventing cellulitis in older people with persistent lower limb oedema. *Br J Nurs* 2004; 13: 725–732.
86. Bonnetblanc JM and Bédane C. Erysipelas: recognition and management. *Am J Clin Dermatol* 2003; 4: 157–163.
87. Villefrance M, Høgh A and Kristensen LH. Compression is important in erysipelas treatment [In Danish]. *Ugeskr Laeger* 2017; 179: pii: V04170284.
88. Ligi D, Croce L and Mannello F. Inflammation and compression: state of the art. *Veins and Lymphatics* 2016; 5: 5980.
89. European Commission. Consumer Goods – Pharmaceuticals: a guideline on summary of product characteristics (SmPC): section 4.8: undesirable effects, www.ema.europa.eu/en/documents/presentation/presentation-section-48-undesirable-effects_en.pdf (2009, accessed July 2019).

INFORMATIONEN ZU DEN AUTOREN

- 1 Department of Dermatology, University of Bonn, Bonn, Germany
- 2 Emeritus, Department of Dermatology, Medical University of Vienna, Vienna, Austria
- 3 Center for Vein Restoration, Mesa, AZ, USA
- 4 Department of Surgery, Harborview Medical Center, Seattle, WA, USA
- 5 Angiology Department, Clinica MD, Barbantini, Lucca, Italy
- 6 Josef Pflug Vascular Laboratory, West London Vascular and Interventional Centre, Ealing Hospital & Imperial College, London, UK
- 7 Centre de Recherche Universitaire de La Lechère, Equipe THEMAS, Université Joseph Fourier, Grenoble, France
- 8 SIGVARIS Management AG, Winterthur, Switzerland
- 9 Department of Dermatology, University of Greifswald, Greifswald, Germany
- 10 Department of General Surgery, Vascular Surgery, Angiology and Phlebology, Medical University of Silesia, Katowice, Poland
- 11 Department of Dermatology, University Hospital of Zurich, Zurich, Switzerland
- 12 First Choice Vascular, Swastik Society, Ahmedabad, India
- 13 Department of Podiatric Surgery & Applied Biomechanics, Rosalind, Franklin University of Medicine and Science, North Chicago, IL, USA
- 14 Department of Vascular Surgery, The University of Chicago, Chicago, IL, USA
- 15 Department of Vascular Surgery, University of Michigan, Ann Arbor, MI, USA
- 16 Practice for Internal Medicine and Vascular Diseases, Halle, Germany

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