SECURITE PRO

by SIGVARIS

Gebrauchsanweisung	2
Instructions for use	4
Istruzioni per l'uso	6
Mode d'emploi	8

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SECURITE PRO

Intended Purpose

Support stocking to support venous return. The product is intended for multiple uses and for application on the legs (lower leg).

The Effects of Medical Support Stockings

The noticeable and controlled compression effect supports the return of blood to the heart, which helps to prevent tired and heavy legs.

Target Group

Support stockings are designed for people who have either a standing or sedentary job, who travel regularly (car, train, and plane), and who want to take preventive action to protect their health.

Indications & Contraindications

Indications

Prevention of tired and heavy legs during long periods of sitting or standing, even without the presence of varicose veins, prevention of latent swelling, subsequent treatment, and prevention following successful varicose vein treatment.

Contraindications

Medical support stockings must not be worn in the following cases:

- Advanced peripheral arterial occlusive disease (if one of the following parameters applies: ABPI <0.5, ankle blood pressure <60 mmHg, toe blood pressure <30 mmHg, or TcPO₂ <20 mmHg instep).
- Decompensated heart failure (NYHA III + IV)
- Septic thrombophlebitis
- Phlegmasia cerulea dolens

Relative contraindications for which the benefits and risks must be carefully weighed up:

- Pronounced exuding dermatoses
- Intolerance to compression material
- Severe sensory disturbances of the extremity
- Advanced peripheral neuropathy (e.g., in diabetes mellitus)
- Primary chronic polyarthritis
- Arterial perfusion disorder (ABPI between 0.5 and 0.8)

Risks and Side Effects

Medical support stockings can cause skin necrosis and pressure damage to the peripheral nerves if handled improperly. If problems should arise, the support stockings should no longer be worn.

Itching, scaling, and inflammation may occur under the support stockings in cases of sensitive skin, so adequate skin care under the compression aids is advised.

One of the yarns used is treated with zinc oxide.

Special Instructions

Do not cut loose ends on the insides of the stocking.

Protect the support stockings from ointments, oils, and grease to preserve their service life.

Support stockings must not be worn over open wounds.

Failure to comply with this information may reduce the safety and effectiveness of the product.

All serious incidents that occur when using the product should be reported to the manufacturer and the competent supervisory authority.

Application

Sharp fingernails and rings can cause defects.

Donning

1. The best time to put the stockings on is after waking up in the morning, when the swelling in your legs has gone down. For easier application, make sure that your legs are dry.
2. Slip the stocking evenly over your heel. Do not pull, but spread the fabric evenly over the leg using the gloves or your hand. While doing so, check for a proper fit at the toe and heel.
3. If it needs repositioning, work the fabric upwards towards the knee with the palm of your hand and spread evenly throughout.

Correct Fit

Knee high stockings:

- 2 finger widths below the back of the knee

Removing

1. Take hold of the stockings at the upper edge with gloves or the palm of your hand and pull downwards towards the heel.
2. Work the fabric over the heel with the palm of your hand.

For advice on donning and removing our support stockings, we recommend watching the videos available at www.sigvaris.com/videos.

**Storage, Wearing Period, and Service Life**

We recommend wearing the support stockings daily from morning to night. Wearing and washing your support stocking every day may cause the medical pressure gradient and elasticity of your support stocking to diminish.

The shelf life date is indicated on the packaging with the hourglass symbol. We guarantee maximum effect for the subsequent six months of wearing.

Care and Storage

Follow the instructions for care on the packaging and/or the textile labels.

Wash inside a mesh laundry bag. The stocking must be stored at room temperature and in a dry place away from light. Keep away from direct heat.

Correct Care**Material Abbreviations:**

EL	elastane
PA	polyamide
AR	aramid

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SIGVARIS and SIGVARIS GROUP are registered trademarks of SIGVARIS AG,
© 2022 Copyright by SIGVARIS AG, CH-9014 St. Gallen/Switzerland
SIGVARIS S.r.l., Via Galileo Galilei 2 / 2, 31057 Silea (TV), **Italy**
SIGVARIS S.A., ul. Powidzka 50, Gutowo Małe, 62-300, Września, **Poland**
SIGVARIS (Australia) Pty. Ltd., 8/53 Stanley Street, West Melbourne VIC 3003, **Australia**
SIGVARIS (SHANGHAI) TRADING CO LTD, Pudong New District, Unit 207, Building 1,