



## COMFORT 850

for women & men

### Compression Level

20-30 mmHg

### Features and Benefits



Cushioned  
Sole



Durable



Moisture  
Wicking

- NILIT Aquarius fibers provide moisture wicking to keep feet dry and comfortable
- Cushioned soles to absorb impact and help prevent blisters
- Smooth, flat toe seam
- Suitable for uniforms requiring solid black color



### Styles and Colours

Closed Toe = ●

BLACK  
99

Calf (C)	●
----------	---

Graduated Compression  
for Men & Women



sigvaris

MOTION



**BLACK**  
**99**

**COMFORT 850**

for women & men

**Compression Level**

20-30 mmHg