

Patient Name:	Contact Name:
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:

Product Information

Product includes one Chipsleeve Standard Arm, one pair of Cotton Liners, and one black Oversleeve.

<input type="checkbox"/> Right Arm	Size:	Item #:
Colour: Black	Length:	Quantity:
<input type="checkbox"/> Left Arm	Size:	Item #:
Colour: Black	Length:	Quantity:

Oversleeve Arm Size Chart

Colour	Small	Medium	Large	X Large	XX Large
Black	1801-ARR	1802-ARR	1803-ARR	1804-ARR	1805-ARR

Circumference

Left Right

C _____

B1 _____

B _____


A _____

D _____

D Length

Short: 38–42.9cm
 Regular: 43–47.9cm
 Long: 48–53cm

Measure length along dorsal aspect



Chipsleeve Standard Arm Size Chart (A–C)

(Available in Black Only)

	Small	Medium	Large	X Large	XX Large
C	23–32cm	28–37cm	33–43cm	39–49cm	44–54cm
B1	22–30cm	26–34cm	30–39cm	35–44cm	39–49cm
B	20–27cm	24–31cm	28–35cm	32–39cm	36–43cm
A	14–17cm	15–18cm	16–19cm	18–20cm	19–22cm

Left

Short	2631-ARS-L	2632-ARS-L	2633-ARS-L	2634-ARS-L	2635-ARS-L
Regular	2631-ARR-L	2632-ARR-L	2633-ARR-L	2634-ARR-L	2635-ARR-L
Long	2631-ARL-L	2632-ARL-L	2633-ARL-L	2634-ARL-L	2635-ARL-L

Right

Short	2631-ARS-R	2632-ARS-R	2633-ARS-R	2634-ARS-R	2635-ARS-R
Regular	2631-ARR-R	2632-ARR-R	2633-ARR-R	2634-ARR-R	2635-ARR-R
Long	2631-ARL-R	2632-ARL-R	2633-ARL-R	2634-ARL-R	2635-ARL-R