

| | |
|----------------------|---------------|
| Patient Name: | Contact Name: |
| Account Name: | Account #: |
| Account Phone #: | P.O. #: |
| Ship Name & Address: | Date: |

Product Information

Product includes one Comprefit Standard Calf wrap, one flat foam liner, one Compreboot® Standard Foot wrap, and 1 pair of Basic Liners. The Compreboot Standard Foot wrap included matches in size of the Comprefit Standard Calf wrap selected. **Note:** If Compreboot needed is a different size, select the Comprefit Plus Calf with separate foot piece.

| | | |
|-------------------------------------|---------|-----------|
| <input type="checkbox"/> Right Calf | Size: | Item #: |
| Colour: Black | Length: | Quantity: |
| <input type="checkbox"/> Left Calf | Size: | Item #: |
| Colour: Black | Length: | Quantity: |

Companion Products (sold separately)

| | |
|---|-----------|
| <input type="checkbox"/> Comprefit Strap Extender | Item #: |
| Colour: <input type="checkbox"/> Black | Quantity: |

Note: Liner and additional Compreboot options available, each sold separately.

Compreboot Standard Foot Size Chart (I, J, K)

| | Small | | Medium / Large | | X Large / XX Large | |
|--------------|----------|----------|----------------|----------|--------------------|----------|
| | Regular | Long | Regular | Long | Regular | Long |
| I | 28–36cm | 28–36cm | ≤39cm | ≤39cm | ≤44cm | ≤44cm |
| J | 22–26cm | 22–26cm | ≤30cm | ≤30cm | ≤33cm | ≤33cm |
| K | 14–18cm | 19–23cm | 18–20cm | 23–25cm | 20–22cm | 25–27cm |
| Black | 1801-BTR | 1801-BTL | 1803-BTR | 1803-BTL | 1805-BTR | 1805-BTL |

Comprefit Strap Extender Chart

| | 1 Strap | 10 Straps |
|--------------|----------|-----------|
| Black | 1100-SE1 | 1100-SEB |

Circumference


Left Right

| | | |
|----|-------|-------|
| C1 | _____ | _____ |
| C | _____ | _____ |
| B | _____ | _____ |
| A | _____ | _____ |
| I | _____ | _____ |
| J | _____ | _____ |

G Length

Regular: 30–35.9cm
Tall: 36–42cm

| | | |
|---|-------|-------|
| G | _____ | _____ |
| K | _____ | _____ |



Measuring Instructions

Foot Size:

- I:** Measure circumference of ankle bend and heel
- J:** Measure circumference across metatarsal heads
- K:** Measure medial length from heel to 1st metatarsal head

Comprefit Standard Calf Regular Size Chart (A–C1)

| Regular | | Small | Medium | Large | X Large | XX Large |
|--------------|----------|----------|----------|----------|----------|----------|
| | C1 | 24–34cm | 29–39cm | 34–44cm | 43–53cm | 50–60cm |
| C | 24–34cm | 29–39cm | 34–44cm | 43–53cm | 50–60cm | |
| B | 19–29cm | 24–34cm | 28–38cm | 36–46cm | 39–50cm | |
| A | 15–24cm | 16–25cm | 20–31cm | 27–37cm | 28–38cm | |
| Black | 1101-BPR | 1102-BPR | 1103-BPR | 1104-BPR | 1105-BPR | |

Comprefit Standard Calf Tall Size Chart (A–C1)

| Tall | | Small Tall | Medium Tall | Large Tall | X Large Tall | XX Large Tall |
|--------------|----------|------------|-------------|------------|--------------|---------------|
| | C1 | 24–34cm | 29–39cm | 34–44cm | 43–53cm | 50–60cm |
| C | 21–31cm | 26–36cm | 30–40cm | 39–49cm | 45–55cm | |
| B | 16–26cm | 20–30cm | 25–35cm | 31–41cm | 35–45cm | |
| A | 15–24cm | 16–25cm | 20–31cm | 27–37cm | 28–38cm | |
| Black | 1101-BPT | 1102-BPT | 1103-BPT | 1104-BPT | 1105-BPT | |