

Patient Name:	Contact Name:
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:

Product Information

Product includes one Comprefit Standard Calf wrap, one Compreboot® Standard Foot wrap, and 1 pair of Basic Liners. The Compreboot Standard Foot wrap included matches the size of the Comprefit Standard Calf wrap selected. **Note:** If Compreboot needed is a different size, select the Comprefit Standard Calf with separate foot piece.

<input type="checkbox"/> Right Calf	Size:	Item #:
Colour: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length:	Quantity:
<input type="checkbox"/> Left Calf	Size:	Item #:
Colour: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length:	Quantity:
<input type="checkbox"/> Comprefit Strap Extender	Item #:	
Colour: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Quantity:	

Note: Liner and additional Compreboot options available, each sold separately.

Compreboot Standard Foot Size Chart (I, J, K)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28–36cm	28–36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22–26cm	22–26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14–18cm	19–23cm	18–20cm	23–25cm	20–22cm	25–27cm
Black	1801-BTR	1801-BTL	1803-BTR	1803-BTL	1805-BTR	1805-BTL
Beige	1811-BTR	1811-BTL	1813-BTR	1813-BTL	1815-BTR	1815-BTL

Comprefit Strap Extender Chart

	1 Strap	10 Straps
Black	1100-SE1	1100-SEB
Beige	1110-SE1	1110-SEB

Circumference


Left Right

C1	_____	_____
C	_____	_____
B	_____	_____
A	_____	_____
I	_____	_____
J	_____	_____

G Length

Regular: 30–35.9cm
Tall: 36–42cm

G	_____	_____
K	_____	_____



Measuring Instructions

Foot Size:

- I:** Measure circumference of ankle bend and heel
- J:** Measure circumference across metatarsal heads
- K:** Measure medial length from heel to 1st metatarsal head

Comprefit Standard Calf Regular Size Chart (A–C1)

	Small	Medium	Large	X Large	XX Large
	Regular				
C1	29–39cm	34–44cm	39–49cm	48–58cm	55–65cm
C	29–39cm	34–44cm	39–49cm	48–58cm	55–65cm
B	24–34cm	29–39cm	33–43cm	41–51cm	44–55cm
A	20–29cm	21–30cm	25–36cm	32–42cm	33–43cm
Black	1101-BKR	1102-BKR	1103-BKR	1104-BKR	1105-BKR
Beige	1111-BKR	1112-BKR	1113-BKR	1114-BKR	1115-BKR

Comprefit Standard Calf Tall Size Chart (A–C1)

	Small Tall	Medium Tall	Large Tall	X Large Tall	XX Large Tall
	Tall				
C1	29–39cm	34–44cm	39–49cm	48–58cm	55–65cm
C	26–36cm	31–41cm	35–45cm	44–54cm	50–60cm
B	21–31cm	25–35cm	30–40cm	36–46cm	40–50cm
A	20–29cm	21–30cm	25–36cm	32–42cm	33–43cm
Black	1101-BKT	1102-BKT	1103-BKT	1104-BKT	1105-BKT
Beige	1111-BKT	1112-BKT	1113-BKT	1114-BKT	1115-BKT