

Patient Name:	Contact Name:
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:

Product Information

Product includes one Compreflex Reduce Calf and one pair of Cotton Liners.

<input type="checkbox"/> Right Calf	Size:	Item #:
Colour: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length:	Quantity:
<input type="checkbox"/> Left Calf	Size:	Item #:
Colour: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length:	Quantity:

Companion Products (sold separately)

<input type="checkbox"/> Compreflex Reduce Foot	Colour: <input type="checkbox"/> Black <input type="checkbox"/> Beige
Size:	Item #: Quantity:
Size:	Item #: Quantity:

Note: Compreboot® options available, each sold separately.

Measuring Instructions

Compreflex Reduce Calf:


Measure circumferences where indicated, beginning 5cm above the lateral malleolus (outer ankle bone), continuing to 25cm. Record in the corresponding lines A and C (Left/Right). Measure length from the lateral malleolus (outer ankle bone) to lateral aspect of the fibular head, and record in box G.

Circumference

	Left	Right
C	_____	_____
A	_____	_____
I	_____	_____
J	_____	_____

G Length
Regular: 30–33.9cm
Tall: 34–42cm

G	_____	_____
K	_____	_____



Measuring Instructions

Foot Size:

- I:** Measure circumference of ankle bend and heel
- J:** Measure circumference across metatarsal heads
- K:** Measure medial length from heel to 1st metatarsal head

Compreflex Reduce Calf Size Chart (A–C)

	Range 1	Range 3
C	<65cm	<90cm
A	<55cm	<75cm

Black		
Regular	2401-BKR	2403-BKR
Tall	2401-BKT	2403-BKT
Beige		
Regular	2411-BKR	2413-BKR
Tall	2411-BKT	2413-BKT

Compreflex Reduce Foot Size Chart (I, J, K)

	Range 1	Range 3
I	30–44cm	39–54cm
J	24–33cm	30–43cm
K	16–25cm	20–27cm
Black		
	2801-BTL	2803-BTL
Beige		
	2811-BTL	2813-BTL