

Patient Name:	Contact Name:
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:

### Circumference


Left    Right

I    \_\_\_\_\_

J    \_\_\_\_\_

### Length

K    \_\_\_\_\_



### Measuring Instructions

#### Foot Size:

- I:** Measure circumference of ankle bend and heel
- J:** Measure circumference across metatarsal heads
- K:** Measure medial length from heel to 1st metatarsal head

### Product Information

Product includes one Compreflex Reduce Foot.

<input type="checkbox"/> Right Foot	Size:	Item #:
Colour: <input type="checkbox"/> Black <input type="checkbox"/> Beige		Quantity:
<input type="checkbox"/> Left Foot	Size:	Item #:
Colour: <input type="checkbox"/> Black <input type="checkbox"/> Beige		Quantity:

**Note:** Compreflex Reduce Calf options available, each sold separately.

### Compreflex Reduce Foot Size Chart (I, J, K)

	Range 1	Range 3
I	30-44cm	39-54cm
J	24-33cm	30-43cm
K	16-25cm	20-27cm
<b>Black</b>	2801-BTL	2803-BTL
<b>Beige</b>	2811-BTL	2813-BTL