

Patient Name:	Contact Name:
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:

Circumference


Left Right

I _____

J _____

Length

K _____



Measuring Instructions

Foot Size:

- I:** Measure circumference of ankle bend and heel
- J:** Measure circumference across metatarsal heads
- K:** Measure medial length from heel to 1st metatarsal head

Product Information

Compreboot® Lite Foot Right Left

Colour: Black Size: Item #: Quantity:

Compreboot Standard Foot Right Left

Colour: Black Beige Item #: Size: Quantity:

Compreboot Plus Foot Right Left

Colour: Black Beige Item #: Size: Quantity:

Coolflex Standard Foot Right Left

Colour: Black Size: Item #: Quantity:

Compreboot Lite Foot Size Chart (I, J, K)

(Available in Black Only)

	Small / Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long
I	28-36cm	28-36cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤33cm	≤33cm
K	14-18cm	19-23cm	20-22cm	25-27cm
Black	1803-BLR	1803-BLL	1805-BLR	1805-BLL

Coolflex Standard Foot Size Chart (I, J, K)

(Available in Black Only)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28-39cm	N/A	30-33cm	N/A	34-49cm	N/A
J	22-26cm	N/A	25-30cm	N/A	28-33cm	N/A
K	14-18cm	19-21cm	18-20cm	20-22cm	20-22cm	23-25cm
Black						
Right	2101-BTR-R	2101-BTL-R	2103-BTR-R	2103-BTL-R	2105-BTR-R	2105-BTL-R
Left	2101-BTR-L	2101-BTL-L	2103-BTR-L	2103-BTL-L	2105-BTR-L	2105-BTL-L

Compreboot Standard Foot Size Chart (I, J, K)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28-36cm	28-36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm
Black	1801-BTR	1801-BTL	1803-BTR	1803-BTL	1805-BTR	1805-BTL
Beige	1811-BTR	1811-BTL	1813-BTR	1813-BTL	1815-BTR	1815-BTL

Compreboot Plus Foot Size Chart (I, J, K)

(Available in Black Only)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28-36cm	28-34cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-25cm	≤30cm	≤28cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm
Black	1801-BPR	1801-BPL	1803-BPR	1803-BPL	1805-BPR	1805-BPL