

Patient Name:	Contact Name:
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:

### Product Information

Product includes one Medaform Standard Calf, one Compreboot® Plus Foot wrap, and one pair of Basic Liners. Indicate below the required Compreboot Plus Foot size.

<input type="checkbox"/> Right Calf	Size:	Item #:
Colour: Blue & Black	Length:	Quantity:
<input type="checkbox"/> Right Foot	Size:	Item #:
Colour: Black	Length:	Quantity:
<input type="checkbox"/> Left Calf	Size:	Item #:
Colour: Blue & Black	Length:	Quantity:
<input type="checkbox"/> Left Foot	Size:	Item #:
Colour: Black	Length:	Quantity:

**Note:** Liner and additional Compreboot options available, each sold separately.

### Compreboot Plus Foot Size Chart (I, J, K) (Available in Black Only)


	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28-36cm	28-34cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-25cm	≤30cm	≤28cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm
<b>Black</b>	1801-BPR	1801-BPL	1803-BPR	1803-BPL	1805-BPR	1805-BPL

### Circumference

	Left	Right
C1	_____	_____
C	_____	_____
B	_____	_____
A	_____	_____
I	_____	_____
J	_____	_____

**G Length**  
 Short: 25-29.9cm  
 Regular: 30-35.9cm  
 Tall: 36-42cm

G	_____
K	_____



### Measuring Instructions

#### Foot Size:

- I:** Measure circumference of ankle bend and heel
- J:** Measure circumference across metatarsal heads
- K:** Measure medial length from heel to 1st metatarsal head

### Medaform Standard Calf Size Chart (A-C1) (Available in Blue & Black Only)

	Small	Medium	Large	X Large	XX Large
C1	32-42cm	38-48cm	42-52cm	51-61cm	58-68cm
C	29-39cm	34-44cm	39-49cm	48-58cm	55-65cm
B	24-34cm	29-39cm	33-43cm	41-51cm	44-55cm
A	20-29cm	21-30cm	25-36cm	32-42cm	33-43cm

		Blue & Black				
Left	Short	1201-BKS-L	1202-BKS-L	1203-BKS-L	1204-BKS-L	1205-BKS-L
	Regular	1201-BKR-L	1202-BKR-L	1203-BKR-L	1204-BKR-L	1205-BKR-L
	Tall	1201-BKT-L	1202-BKT-L	1203-BKT-L	1204-BKT-L	1205-BKT-L
Right	Short	1201-BKS-R	1202-BKS-R	1203-BKS-R	1204-BKS-R	1205-BKS-R
	Regular	1201-BKR-R	1202-BKR-R	1203-BKR-R	1204-BKR-R	1205-BKR-R
	Tall	1201-BKT-R	1202-BKT-R	1203-BKT-R	1204-BKT-R	1205-BKT-R