

Patient Name:		Contact Name:
Account Name:		Account #:
Account Phone #:		P.O. #:
Item #: 2739-AR	Height:	Weight:
Ship Name & Address:		Date:

Product Information

Product includes one Chipsleeve Custom Arm, one pair of Cotton Liners, and one black Oversleeve.

- | | |
|------------------------------------|---------------|
| <input type="checkbox"/> Right Arm | Colour: Black |
| <input type="checkbox"/> Left Arm | Colour: Black |

Important

Measuring for this garment requires special measurements critical to proper fit. Before you begin, please contact Customer Care, or your Territory Manager for more information.

Supplies Needed

- Cell phone with camera. Photos of the arm with measurement markings must be emailed to: Canadaorders@sigvaris.com
- Measuring instructions and forms.
- SIGVARIS GROUP Measuring tape and body pen (or eyeliner pencil).
- Signed Custom Order Terms & Conditions Form (include with your order).

Circumference

Left Right

C _____

B1 _____

B _____

A1 _____


A _____

E _____

D Length

Measure length along dorsal aspect

D _____



MEASURING INSTRUCTIONS

With patient seated, place the arm extended and elbow bent slightly, with the palm down, on a flat surface.

CIRCUMFERENCES

- Measure circumference at palm and record on line E.
- Measure circumference at wrist, and record on line A.
- Mark dorsal aspect at distal edge of tape.
- Measure circumference at forearm, and record on line A1.
- Measure circumference at elbow, and record on line B.
- Measure circumference at mid bicep, and record on line B1.
- Measure circumference at axilla, and record on line C.
- Mark dorsal aspect at proximal edge of tape.

LENGTH

Measure length at dorsal aspect from mark at Point A to mark at Point C and record in box D.