


Patient Name:	Contact Name:
Account Name:	Account #:
Account Phone #:	P.O. #:
Height:	Weight:
Ship Name & Address:	Date:



Product Information

Product includes one Chipvest Custom.

Custom Chipvest Item #: **2239-VS**

- Full (Bilateral)
- Right Side (Unilateral)
- Left Side (Unilateral)

Measuring Instructions

Step One:

Measure the length in centimeters from the suprasternal notch to the waist.

Step Two:

Measure circumferences in centimeters.

Custom Size Measurements

Length	Centimeters
L	
Circumference	Centimeters
Axilla	
Largest Chest	
Xyphoid Process	
Waist	
Mid-Hip	

Caution: Drawstring contains natural rubber latex.

Important

Measuring for this garment requires special measurements critical to proper fit. Before you begin, please contact Customer Care, or your Territory Manager for more information.

Supplies Needed

- Cell phone with camera. Photos of the torso with measurement markings must be emailed to: **Canadaorders@sigvaris.com**
- Measuring instructions and forms.
- SIGVARIS GROUP Measuring tape and body pen (or eyeliner pencil).
- Signed Custom Order Terms & Conditions Form (include with your order).