

Patient Name:	Contact Name:
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:

Product Information

Product includes one Medaform Full Leg, one Compreboot® Plus Foot, and one pair of Cotton Liners. Indicate below the required Compreboot Plus Foot size.

<input type="checkbox"/> Right Leg	Size:	Item #:
Colour: Blue & Black	Length:	Quantity:
<input type="checkbox"/> Right Foot	Size:	Item #:
Colour: Black	Length:	Quantity:
<input type="checkbox"/> Left Leg	Size:	Item #:
Colour: Blue & Black	Length:	Quantity:
<input type="checkbox"/> Left Foot	Size:	Item #:
Colour: Black	Length:	Quantity:

Note: Liner and additional Compreboot options available, each sold separately.

Compreboot Plus Foot Size Chart (I, J, K) (Available in Black Only)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28-36cm	28-34cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-25cm	≤30cm	≤28cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm
Black	1801-BPR	1801-BPL	1803-BPR	1803-BPL	1805-BPR	1805-BPL


Circumference

	Left	Right
F	_____	_____
E	_____	_____
D	_____	_____
Ø Point	_____	_____
C1	_____	_____
C	_____	_____
B	_____	_____
A	_____	_____
I	_____	_____
J	_____	_____

H Length

Short: 61-70.9cm
Regular: 71-80.9cm
Tall: 81-91cm

H	_____	_____
K	_____	_____



Measuring Instructions

Foot Size:

- I:** Measure circumference of ankle bend and heel
- J:** Measure circumference across metatarsal heads
- K:** Measure medial length from heel to 1st metatarsal head

Medaform Full Leg Size Chart (A-F) (Available in Blue & Black Only)

	Small	Small Max	Medium	Medium Max	Large	Large Max	X Large
F	48-58cm	56-66cm	56-66cm	64-74cm	64-74cm	74-84cm	74-84cm
E	43-53cm	51-61cm	51-61cm	58-68cm	58-68cm	68-78cm	68-78cm
D	38-48cm	46-56cm	46-56cm	53-63cm	53-63cm	63-73cm	63-73cm
C1	32-42cm		37-47cm		42-52cm		51-61cm
C	29-39cm		34-44cm		39-49cm		48-58cm
B	24-34cm		29-39cm		33-43cm		41-51cm
A	20-29cm		21-30cm		25-36cm		32-42cm

Blue & Black

		Blue & Black								
		Short	Regular	Tall	Short	Regular	Tall	Short	Regular	Tall
Left	Short	1201-THS-L	1201-TMS-L	1202-THS-L	1202-TMS-L	1203-THS-L	1203-TMS-L	1204-THS-L		
	Regular	1201-THR-L	1201-TMR-L	1202-THR-L	1202-TMR-L	1203-THR-L	1203-TMR-L	1204-THR-L		
	Tall	1201-THT-L	1201-TMT-L	1202-THT-L	1202-TMT-L	1203-THT-L	1203-TMT-L	1204-THT-L		
Right	Short	1201-THS-R	1201-TMS-R	1202-THS-R	1202-TMS-R	1203-THS-R	1203-TMS-R	1204-THS-R		
	Regular	1201-THR-R	1201-TMR-R	1202-THR-R	1202-TMR-R	1203-THR-R	1203-TMR-R	1204-THR-R		
	Tall	1201-THT-R	1201-TMT-R	1202-THT-R	1202-TMT-R	1203-THT-R	1203-TMT-R	1204-THT-R		