


Patient Name:	Contact Name:
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:

### Circumference

Head \_\_\_\_\_

Neck \_\_\_\_\_



### Product Information

Product includes one Neck and Mandible Pad.

Neck & Mandible      Colour: Black

Size:      Quantity:      Item #:

### Measuring Instructions

#### Head:

Measure the circumference by encircling the head; beginning at the top, moving below the chin; and finishing at the top of the head.

#### Neck:

Measure circumference at largest part of neck.

### Neck & Mandible Pad Sizing Chart (Available in Black Only)

	Small	Medium	Large	X Large
Head	46–56cm	52–62cm	56–66cm	58–68cm
Neck	28–33cm	33–38cm	38–43cm	43–58cm
<b>Black</b>	2231-NM	2232-NM	2233-NM	2234-NM