

Patient Name:	Contact Name:
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:

Product Information

Product includes one Compreflex Reduce Knee wrap and one Butterfly strap

<input type="checkbox"/> Right Knee	Size:	Item #:
Colour: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length:	Quantity:
<input type="checkbox"/> Left Knee	Size:	Item #:
Colour: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length:	Quantity:

Companion Products (sold separately)

<input type="checkbox"/> Compreflex Reduce Thigh	Colour: <input type="checkbox"/> Black <input type="checkbox"/> Beige
Size:	Item #: Quantity:
Size:	Item #: Quantity:

Measuring Instructions

Compreflex Reduce Knee:

Measure circumferences where indicated, beginning 5cm below the lateral aspect fibular head. Continue 5cm above the lateral aspect top of patella. Record in the corresponding lines C1 and D.

Compreflex Reduce Thigh:

Measure circumferences where indicated on the measurement chart, beginning 5cm above the top of patella and continue 23cm up the thigh, just below the gluteal fold. Record in the corresponding lines D and F. Measure length from the lateral aspect top of patella to the gluteal fold, and record in box G1.

Circumference


	Left	Right
F	_____	_____
E	_____	_____
D	_____	_____
Ø Point	_____	_____
C1	_____	_____

G1 Length

Regular: 20–24.9cm
Tall: 25–30cm

Top of patella to gluteal fold

C1 _____



Compreflex Reduce Knee Size Chart

	Range 1	Range 3
D	<75cm	<115cm
C1	<65cm	<90cm
Black	2401-KP	2403-KP
Beige	2411-KP	2413-KP

Compreflex Reduce Thigh Size Chart

		Range 1	Range 3
	F	<80cm	<125cm
	D	<65cm	<110cm
	Black		
Left	Regular	2401-TCR-L	2403-TCR-L
	Tall	2401-TCT-L	2403-TCT-L
Right	Regular	2401-TCR-R	2403-TCR-R
	Tall	2401-TCT-R	2403-TCT-R
	Beige		
Left	Regular	2411-TCR-L	2413-TCR-L
	Tall	2411-TCT-L	2413-TCT-L
Right	Regular	2411-TCR-R	2413-TCR-R
	Tall	2411-TCT-R	2413-TCT-R