

COMPREFLEX® Contain

If you experience any difficulty finding this item on your prescribing system, please contact Customer Services on the phone number at the bottom of the page

Patient and referrer information

Patient Name / NHS number		Patient's DOB / /	
Clinician name		Date / /	
Email	Tel		

Sizing chart (circumference measurements are shown in cm)

COMPREFLEX® Contain Thigh

	S	M	L	XL	XXL
Above knee (E)	38-48	46-56	53-63	61-71	69-79
Thigh (F)	43-53	51-61	58-68	66-76	74-84
Below gluteal fold (G)	48-58	56-66	64-74	72-82	80-90

Length

	Regular	Tall
Thigh (L2)	20-24.9	25-30

COMPREFLEX® Contain Knee

	S	M	L	XL	XXL
Above Calf (D)	29-39	34-44	39-49	46-56	54-64
Above knee (E)	38-48	46-56	53-63	61-71	69-79

COMPREFLEX® Contain Calf

	S	M	L	XL	XXL
Ankle (B)	19-26	21-30	26-36	31-41	36-46
Calf (C)	29-39	34-44	40-50	46-56	54-64

Length

	Regular	Tall
Calf (L1)	30-33.9	34-42

Thigh and knee

	Left	Right
G		
F		
E		
L2		
D		
E		

Calf

	Left	Right
C		
B		
L1		

Product description

Ordering details

COMPREFLEX® Contain Thigh and Knee

	COMPREFLEX® Contain Thigh				COMPREFLEX® Contain Knee
	Regular		Tall		
	Left	Right	Left	Right	
S	☐	☐	☐	☐	☐
M	☐	☐	☐	☐	☐
L	☐	☐	☐	☐	☐
XL	☐	☐	☐	☐	☐
XXL					
Quantity					

COMPREFLEX® Contain Calf

	COMPREFLEX® Contain Calf	
	Regular	Tall
S	☐	☐
M	☐	☐
L	☐	☐
XL	☐	☐
XXL	☐	☐
Quantity		

Includes 1 wrap (20-50 mmHg) and 2 Transition liners (10-20 mmHg in the foot and ankle)

Additional liners

	For below knee			For thigh
	Basic liner (only in grey)	Transition liner	Complete liner	Cotton liner (only in white)
S	☐ One size fits all	☐	☐	☐ One size fits all
M		☐	☐	
L		☐	☐	
XL		☐	☐	
XXL		☐	☐	
Quantity				

10-20 mmHg in the foot and ankle

10-20 mmHg in the calf, foot and ankle

Strap extenders

Quantity

Number required	
-----------------	--

Hook extenders

Quantity

Number required	
-----------------	--

Product description

SIGVARIS Britain Ltd

1 Imperial Court, Magellan Close,
Andover, Hampshire, SP10 5NT

01264 326666
sigvaris.com

SIGVARIS GROUP is a registered trademark of SIGVARIS AG, St. Gallen/Switzerland, in many countries worldwide. © 2026 Copyright by SIGVARIS AG, St. Gallen/Switzerland 2026/200

With over 160
years of expertise
within medical
compression wear,
SIGVARIS GROUP
is dedicated to help
people feel their
best. Every day

SIGVARIS
GROUP