

Essential THERMOREGULATING

If you experience any difficulty finding this item on your prescribing system, please contact Customer Services on the phone number at the bottom of the page


Patient and referrer information

Patient Name / NHS number		Patient's DOB
Clinician name		Date
Email	Tel	

Sizing chart (circumference measurements are shown in cm)

Below knee, Thigh, Thigh w/attachment and Tights (panty)

	XXS		XS		S		M		L		XL	
Ankle (cB)	18-20		20-22		22-24		24-26		26-29		29-32	
	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS	STD	PLUS
Calf (cC)	25-32	29-36	29-36	33-40	33-40	35-42	35-42	38-45	38-45	42-49	42-49	43-51
Thigh (cG) with grip top	42-52	47-56	47-56	52-61	52-61	56-65	56-65	61-71	61-71	66-76	66-76	73-85
Thigh (cG) <small>(For thigh w/attachment and tights only)</small>	42-54	47-59	47-59	52-64	52-64	56-68	56-68	61-74	61-74	66-79	66-79	76-85
Hip (cH) <small>(For thigh w/attachment and tights only)</small>	75-90	80-95	80-95	85-100	85-100	95-110	95-110	105-120	105-120	115-130	115-130	125-140
Waist (cT) <small>(For thigh w/attachment and tights only)</small>	60-75	65-80	65-80	75-90	75-90	85-100	85-100	95-110	95-110	105-120	105-120	115-130



Length	Short	Normal	Long
Below knee (lD)	31-34	35-40	41-45
Thigh (lG)	56-64	65-74	75-84
Tights (lG)	56-64	65-74	75-84

Patient measurements		cB	cC	cG	cH	cT	Length
Left							
Right							

Ordering details

Product style and size

		Below knee			Below knee Maxi		Thigh			Thigh w/attachment						Tights			
		Short	Normal	Long	Normal	Long	Short	Normal	Long	Short		Normal		Long		Short	Normal	Long	
										Left	Right	Left	Right	Left	Right				
XXS	STD																		
	PLUS																		
XS	STD																		
	PLUS																		
S	STD																		
	PLUS																		
M	STD																		
	PLUS																		
L	STD																		
	PLUS																		
XL	STD																		
	PLUS																		
Quantity																			

Colour

Black	
Nature	

Toe

Closed toe	
Open toe	

Compression level

Class 1 (18-21 mmHg)	
Class 2 (23-32 mmHg)	
Class 3 (34-46 mmHg)	

Grip top
(Below knee in Class 2 and Thigh products)

Knobbed grip top	
Sensinnov grip top	

Sensinnov is only available for thigh

PLEASE NOTE: Not all combinations of Essential THERMOREGULATING (depending on your choice of Product style, size, compression class, colour, grip top or toe) are valid. Please check your prescription order using the latest SIGVARIS Healthcare Professional Guide.

SIGVARIS Britain Ltd
 1 Imperial Court, Magellan Close,
 Andover, Hampshire, SP10 5NT
01264 326666
sigvaris.com

SIGVARIS GROUP is a registered trademark of SIGVARIS AG, St. Gallen/Switzerland, in many countries worldwide. © 2024/2025 Copyright by SIGVARIS AG, St. Gallen/Switzerland

With 160 years of expertise within medical compression wear, **SIGVARIS GROUP** is dedicated to help people feel their best. Every day

