

Style PATTERNS

If you experience any difficulty finding this item on your prescribing system, please contact Customer Services on the phone number at the bottom of the page

Patient and referrer information

Patient Name / NHS number		Patient's DOB
Clinician name		Date
Email	Tel	

Sizing chart (circumference measurements are shown in cm)


Below knee

	S	M	L	XL
Ankle (cB)	19-22	22-24	24-26	26-29
Calf (cC)	32-38	35-41	38-44	41-47


Thigh / Tights (panty)

	S	M	L	XL
Ankle (cB)	19-22	22-24	24-26	26-29
Calf (cC)	32-38	35-41	38-44	41-47
Thigh (cG)	45-58	48-61	51-65	54-68
Hip (cH) (For tights only)	≤ 105	≤ 115	≤ 125	≤ 125


Below knee



Thigh



Tights (panty)



Length

	Normal	Long
Below knee (lD)	≤ 40	≥ 40
Thigh (lG)	≤ 75	≥ 75
Tights (lG)	≤ 75	≥ 75

Patient measurements

	cB	cC	cG	cH	Length
Left					
Right					



Need some help?

Scan the QR code to watch our videos on how to measure

Ordering details

Product style and size

	Below knee		Thigh		Tights	
	Normal	Long	Normal	Long	Normal	Long
S						
M						
L						
XL						
Quantity						

Pattern

Chequered	
Mosaic	
Polka Dot	

Mosaic pattern is only available in tights

Compression level

Class 2 (15-20 mmHg)	
----------------------	--

SIGVARIS Britain Ltd
1 Imperial Court, Magellan Close,
Andover, Hampshire, SP10 5NT

01264 326666
sigvaris.com

SIGVARIS GROUP is a registered trademark of SIGVARIS AG, St. Gallen/Switzerland, in many countries worldwide. © 2024/2025 Copyright by SIGVARIS AG, St. Gallen/Switzerland

With 160 years of expertise within medical compression wear, **SIGVARIS GROUP** is dedicated to help people feel their best. Every day

SIGVARIS
GROUP