



COMPRESLEEVE BK

orders.au@sigvaris.com • P: 1300 SIGVARIS • F: 03 9329 3773

P.O. #	Account #	Contact	Date			
Phone	Patient Name			Age	Height	Weight
Company Name						
Dropship Name & Address						
Product Information						
Leg <input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT	Size	Length	Item #	CompreBoot Size		
Leg <input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT	Size	Length	Item #	CompreBoot Size		

CIRCUMFERENCE

LEFT RIGHT

_____ _____ **C1**

_____ _____ **C**

_____ _____ **B**


_____ _____ **A**

_____ _____ **Ø POINT**

_____ _____ **I**

_____ _____ **J**

_____ _____ **K:**



COMPRESLEEVE - BK SIZING CHART (A-C1, W/ FOOT)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	32–42cm	38–48cm	42–52cm	51–61cm	58–68cm
C	29–39cm	34–44cm	39–49cm	48–58cm	55–65cm
B	24–34cm	29–39cm	33–43cm	41–51cm	44–55cm
A	20–29cm	21–30cm	25–36cm	32–42cm	33–43cm
SHORT	1501-BKS	1502-BKS	1503-BKS	1504-BKS	1505-BKS
REGULAR	1501-BKR	1502-BKR	1503-BKR	1504-BKR	1505-BKR
TALL	1501-BKT	1502-BKT	1503-BKT	1504-BKT	1505-BKT

COMPREBOOT SIZING CHART (I-K)

	SMALL		MED / LARGE		XL / XXL	
	REGULAR	LONG	REGULAR	LONG	REGULAR	LONG
I	28–36cm	28–36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22–26cm	22–26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14–18cm	19–23cm	18–20cm	23–25cm	20–22cm	25–27cm
BLACK	1801-BTR	1801-BTL	1803-BTR	1803-BTL	1805-BTR	1805-BTL
BEIGE	1811-BTR	1811-BTL	1813-BTR	1813-BTL	1815-BTR	1815-BTL