




FOOT COMPRESSION

orders.au@sigvaris.com • P: 1300 SIGVARIS • F: 03 9329 3773

P.O. #	Account #	Contact	Date			
Phone	Patient Name			Age	Height	Weight
Company Name						
Dropship Name & Address						
Product Information						
Foot <input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT						
Size	<input type="checkbox"/> COMPREBOOT LITE					
Length	<input type="checkbox"/> COMPREBOOT					
Item #	<input type="checkbox"/> COMPREBOOT PLUS					

MEASURING INSTRUCTIONS

- I:** Measure circumference of ankle bend and heel.
- J:** Measure circumference across metatarsal heads.
- K:** Measure length on the medial side of the foot from the heel to the 1st metatarsal head.



CIRCUMFERENCE

LEFT RIGHT

_____ I

_____ J

_____ K

SIZING CHART FOR ALL COMPREBOOT PRODUCTS (I-K)

	SMALL		MEDIUM / LARGE		XLARGE / XXLARGE	
	REGULAR	LONG	REGULAR	LONG	REGULAR	LONG
I	28-36cm	28-36cm	<39cm	<39cm	<44cm	<44cm
J	22-26cm	22-26cm	<39cm	<39cm	<33cm	<33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm

COMPREBOOT LITE PRODUCT CODES

	SMALL / MEDIUM	LARGE / X LARGE
REGULAR	1803-BLR	1805-BLR
LONG	1803-BLL	1805-BLL

COMPREBOOT PRODUCT CODES

	SMALL	MEDIUM / LARGE	XLARGE / XXLARGE
BLACK			
REGULAR	1801-BTR	1803-BTR	1805-BTR
LONG	1801-BTL	1803-BTL	1805-BTL
BEIGE			
REGULAR	1811-BTR	1813-BTR	1815-BTR
LONG	1811-BTL	1813-BTL	1815-BTL

COMPREBOOT PLUS PRODUCT CODES

	SMALL	MEDIUM / LARGE	XLARGE / XXLARGE
REGULAR	1801-BPR	1803-BPR	1805-BPR
LONG	1801-BPL	1803-BPL	1805-BPL

For information regarding the **MEDABOOT CUSTOM**, please speak with your SIGVARIS Sales Representative or Customer Care: 1300 SIGVARIS